

***ClearBand* 6x DNA Loading Dye Green**

5x1 ml

Cat No: DLD-6xG

Shipping : Ship at ambient temperature.

Storage : Stable for 1 year at 4-8°C or at room temperature.

General Information

ClearBand 6x DNA Loading Dye **Green** is a pre-mixed buffer, which can be used for loading DNA markers and samples into the wells of agarose and non-denaturing polyacrylamide gels.

It contains two tracking dyes, Xylene Cyanol FF and Orange G, for visualization of DNA migration on the gels. Orange G runs around 50 bp fragment in agarose gels, which provides much better tracking of DNA before it runs off the gel. Besides, Orange G does not result in formation of a dark shadow on the gel as Bromophenol Blue does, which runs around 400 bp fragment precisely where PCR fragments are likely to be.

Protocol

1. Add 1 volume of *ClearBand* 6x DNA Loading Dye **Green** to 5 volume of your DNA sample.
2. Load into the wells of your gel.